



Forever Medical

FORERMED®
www.forarmedchina.com

8L 12L 18L 23L Class B Pulsating Vacuum Table Top Steam Sterilizer





Features

- ☆ Quality stainless steel inner tank, sealed cover
- ☆ Microcomputer control
- ☆ Safety interlocking device
- ☆ With BD test, HELIX test, VACUUM test.
- ☆ LCD display and touch buttons.
- ☆ With standard test interface
- ☆ Open-type top water tank
- ☆ European CLASS B standard.
- ☆ With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti dry burning, and safety door lock.
- ☆ Three times pre-vacuum
- ☆ Buzzer reminder after sterilization, automatic stop.
- ☆ Automatic discharge of cold air, automatic steam exhaust after sterilization.
- ☆ Built-in steam generator to provide saturated steam quickly
- *Optional external printer or built-in printer**
- *Optional USB function and U disk**

Model	YJ-8L	YJ-12L	YJ-18L	YJ-23L
Technical parameter				
Volume	8L	12L	18L	23L
Power	1.5kw	1.8kw	2.2kw	2.2kw
Voltage	110V/220V,60HZ/50HZ			
Rated working pressure	0.22Mpa			
Rated working temperature	134℃			
Chamber volume	Ø170×320	Ø200×355	Ø250×370	Ø250×470
Outer dimension	445×580×400	445×580×400	505×620×450	505×710×450
Packing size (L×W× H)	650×500×450	650×500×450	670×560×510	760×560×510
Gross/net weight	45/42KG	48/44KG	62/53KG	65/58KG